



FLOOR STANDING HIGH SPEED REFRIGERATED CENTRIFUGE BKE-3901

HIGH SPEED REFRIGERATED CENTRIFUGE BKE-3901

FLOOR HIGH SPEED LOW TEMPERATURE CENTRIFUGE



The series is stable but not heavy, excellent compatibility, large capacity with light rotor.

Controlled by microcomputer, AC frequency variable motor drive, without toner power pollution, and longer lifetime. Lower noise and less vibration, equipped with multiple types of rotors for users to choose.

Auto recognition of rotor and speed limit control, to make centrifuge safer.

Imported compressor set, refrigeration and heating double temperature control, unique heat radiation system, to ensure that the cooling effect of the instrument to achieve the best.

Able to configure acceleration/deceleration mode, provide configuration of RPM and RCF, and convert RPM/RCF using one key.

Equipped with over speed, overheat, imbalance and door lid auto lock protection

SPECIFICATIONS

Model	BKE-3901
Max. Speed	21000rpm
Max. RCF	47580xg
Max. Capacity	4x750ml
timing range	1min~23h59min
Temp Range(freezing)	-20°C~40°C
Temp Accuracy(freezing)	±2.0°C
RPM accuracy	±20%
power supply	AC220±22V 50Hz 35A
total power	3500W
Noise level	≤65dB(A)
centrifuge chamber diameter	500mm
dimensions (length x width x height)	825x625x910(mm)
packaging dimensions (length x width x height)	930x730x1010(mm)
Net Weight	310Kg
Gross Weight	345Kg
Alt Name	Floor High speed Low temperature centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	Angle Rotor	16x10ml	21000rpm	47580xg
2		6x50ml	20000rpm	43000xg
3		8x50ml	15000rpm	27720xg
4		4x250ml	14000rpm	29440xg
5		8x100ml	13000rpm	19830xg
6		6x250ml	12000rpm	22470xg
7		6x300ml	10000rpm	18500xg
8		6x500ml	9000rpm	14400xg
9	Swing-out Rotor	4x750ml	4000rpm	3580xg

FEATURES

The series is stable but not heavy, excellent compatibility, large capacity with light rotor.

Controlled by microcomputer, AC frequency variable motor drive, without toner power pollution, and longer lifetime. Lower noise and less vibration, equipped with multiple types of rotors for users to choose.

Auto recognition of rotor and speed limit control, to make centrifuge safer.

Imported compressor set, refrigeration and heating double temperature control, unique heat radiation system, to ensure that the cooling effect of the instrument to achieve the best.

Able to configure acceleration/deceleration mode, provide configuration of RPM and RCF, and convert RPM/RCF using one key.

Equipped with over speed, overheat, imbalance and door lid auto lock protection systems to ensure safety of machine and users

Compact size, equipped with large torque variable-frequency motor to facilitate acceleration and deceleration and enhance efficiency. Higher performance to cost ratio.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com