



LABORATORY VERTICAL AUTOCLAVE BKK-1002

LABORATORY VERTICAL AUTOCLAVE BKK-1002

VERTICAL PRESSURE STEAM STERILIZER



Fully stainless steel structure.
Auto discharge the cool air, and steam discharge automatically after sterilization.
Automatically shut off with beep reminding after sterilization.
Hand wheel type of quick-open door structure.
Over temperature & over pressure auto-protection.
Easy to operate, safe and reliable.
Door safety lock system.
With two stainless steel sterilizing baskets.
LCD display of working status, touch type key.
Safe protection of water lacking.
Optional drying system.

SPECIFICATIONS

Model	BKK-1002
Capacity	50L
Sterilization Temperature	134°C
Working Pressure	0.23 Mpa
Features	
Chamber volume	50L Φ 340x550mm
Max. working pressure	0.23Mpa
Heat average	$\leq \pm 1^{\circ}\text{C}$
Timer range	0-99min or 0-99hour59min
Temperature adjust range	105-134°C
Power	2X1.5KW/AC220V 50HZ
Overall dimension(mm)	520x490x1150
Transport dimension(mm)	580x580x1250
G.W/ N.W	94Kg/74Kg
Alt Name	Vertical Pressure Steam Sterilizer

FEATURES

Fully stainless steel structure.
Auto discharge the cool air, and steam discharge automatically after sterilization.
Automatically shut off with beep reminding after sterilization.
Hand wheel type of quick-open door structure.
Over temperature & over pressure auto-protection.
Easy to operate, safe and reliable.
Door safety lock system.
With two stainless steel sterilizing baskets.
LCD display of working status, touch type key.
Safe protection of water lacking.
Optional drying system.
Self-inflating type seal.

APPLICATIONS

The vertical pressure steam sterilizers are assembled with a heating system, micro-computer controlled system, over heat and over pressure protection system, which are reliable for the sterilizing effects, convenient for operation and energy conservation. They are ideal equipment for clinics, scientific research institutes and other organizations to sterilize surgical instruments, fabrics, glasses, culture media etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com