



## ULTRASONIC HOMOGENIZER BLIH-303

# ULTRASONIC HOMOGENIZER BLIH-303

## ULTRASONIC HOMOGENIZER

Ultrasonic homogenizer series is designed for applying a cavitation effect in a liquid by means of ultrasonic waves. It can be used to disintegrate cells of various animals and plants, viruses, bacteria, tissues as well as to reshape inorganic substances. It can also be used in emulsification, separation, distribution, collection, cleaning and nano-material preparation, dispersion acceleration of chemical reactions, etc.



7 inch large screen LCD touch screen display.  
Microcomputer control, which can store 20 sets of work data.  
Ultrasonic time, ultrasonic power can be set.  
Automatic detection of ultrasonic power to prevent the change of ultrasonic power as the  
The temperature of the sample is changed.  
Integrated temperature control to prevent sample overheating.  
Automatic frequency tracking, automatic fault alarm.  
Single ultrasound time: 0.1-99.9S  
Single interval time: 0.1-99.9S  
Total working time: 1S-99H59M59S(touch type to 99hours)  
Mode of work: clearance or continuity

## SPECIFICATIONS

Model	BLIH-303
Ultrasonic Frequency (kHz)	20-25
Maximum Power	0-500 W
Optional Probe (Φ mm)	2, 3, 6, 8, 10
Crusher Capacity	0.1-400 ml
Sample Temperature Monitor	Yes
Alt Name	Ultrasonic Homogenizer

## FEATURES

7 inch large screen LCD touch screen display.  
Microcomputer control, which can store 20 sets of work data.  
Ultrasonic time, ultrasonic power can be set.  
Automatic detection of ultrasonic power to prevent the change of ultrasonic power as the  
The temperature of the sample is changed.  
Integrated temperature control to prevent sample overheating.  
Automatic frequency tracking, automatic fault alarm.  
Single ultrasound time: 0.1-99.9S  
Single interval time: 0.1-99.9S  
Total working time: 1S-99H59M59S(touch type to 99hours)  
Mode of work: clearance or continuity  
frequency: 19.5-20.5KHz  
Can select computer online function

## APPLICATIONS

It has wide application in teaching, researching, and manufacturing in many fields such as biology, microbiology, physics, zoology, agriculture and pharmacology and has been used by some labs.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)