



## VERTICAL LAMINAR AIRFLOW BLVR-101

# VERTICAL LAMINAR AIRFLOW BLVR-101

## VERTICAL LAMINAR FLOW CABINET



LCD Display, Airflow velocity, UV timer, UV work time, system work time, real time.  
Waterproof Socket, 2 waterproof sockets are located inside the panel, for optimum convenience of using small devices inside the cabinet.  
Side Glass Window, The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.  
Gas Tap  
UV Lamp, Emission of 253.7 nanometers for most efficient decontamination.

## SPECIFICATIONS

Model	BLVR-101
External Size (W*D*H)	1300*750*2040mm
Internal Size (W*D*H)	1200*645*610mm
Work Surface Height	750mm
Display	LCD Display
Aiflow Velocity	Average of 0.3~0.5m/s
Material	Main Body: Cold-rolled steel with anti-bacteria powder coating Work Table: 304 stainless steel Side and Front Windows: 5mm toughened glass anti-UV
Pre-Filter	Polyester fiber, washable
HEPA Filter	99.995% efficiency at 0.3µm
Noise	≤65dB
Front Window	Motorized
Max Opening	430mm
LED Lamp	12W*1
UV Lamp	30W*1 Emission of 253.7 nanometers
Consumption	400W
Waterproof Socket	Two, total loads 500W
Caster	Universal caster with leveling feet
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz
Standard Accessory	LED lamp, UV lamp*2, Base stand, Gas Tap, Waterproof socket*2
Optional Accessory	Electric height adjustable base stand
Gross Weight	228kg
Package Size (W*D*H)	1460*1070*1650mm
Alt Name	Vertical Laminar Flow Cabinet

## FEATURES

LCD Display, Airflow velocity, UV timer, UV work time, system work time, real time.



Waterproof Socket, 2 waterproof sockets are located inside the panel, for optimum convenience of using small devices inside the cabinet.




Side Glass Window, The transparent side glass windows maximize light and visibility inside the cabinet, providing a bright and open working environment.



Gas Tap



UV Lamp, Emission of 253.7 nanometers for most efficient decontamination.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)