



ELECTRIC OPERATION TABLE BOPT-105

ELECTRIC OPERATION TABLE BOPT-105

ELECTRIC OR TABLE

Biolab's electric operation table are cost effective, multi performance mobile operating table that meets the requirement of different surgeons. Easy-to-use electronic controls. Operation table equipped with a table top that rotates 180 degrees, and the chest plate, back plate, head plate, and leg plate on this table can all be easily adjusted, allowing for the ideal positioning of patients undergoing medical imaging tests.



The tabletop adopts high quality of carbon plastic plate, and it can slide, which is convenient for C arm exam and X-ray photography.
 AC and DC, dual remote control to control various actions
 The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam
 Imported brand LINAK motor, with no noise, long life
 Quality stainless steel 304 is used for table support, base plate, column shield and rails on both sides of the table
 Foot pedal brake device which is easy for moving the operating table
 The head board and leg plate are controlled by imported gas spring
 Leg plate is split and detachable
 It has the function of auto-restoration

SPECIFICATIONS

Model	BOPT-105
Table Length	2000±100 mm
Table Width	520±20 mm
Height Adjustment	(700~1150) ±50 mm
Tilt (left)/Tilt (Right)	≥15°/≥15°
Trendelenburg/Rev.-Trend	≥20°/≥20°
Head board upwards	≥50°
Head board downwards	≥90°
Backboard upwards	≥70°
Backboard downwards	≥40°
Waist board lifting	100±20 mm
Leg board upwards	≥15°
Leg board downwards	≥90°
Leg board outwards	≥180°
Lumbar upper corner "∧"	≤150°
Lumbar lower corner "∨"	≤110°
Horizontal and longitudinal movement	400±30 mm
Standard Accessories	A Pair of Shoulder/Waist Support, A pair of Arm Support, A pair of leg Plate, A screen
Main Voltage AC	220V±22V, 50Hz±1Hz
Input Power VA	400±10%
Alt Name	Electric OR table

FEATURES

The tabletop adopts high quality of carbon plastic plate, and it can slide, which is convenient for C arm exam and X-ray photography.
 AC and DC, dual remote control to control various actions
 The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam
 Imported brand LINAK motor, with no noise, long life
 Quality stainless steel 304 is used for table support, base plate, column shield and rails on both sides of the table
 Foot pedal brake device which is easy for moving the operating table
 The head board and leg plate are controlled by imported gas spring

Leg plate is split and detachable
It has the function of auto-restoration
One press to "Λ" or "∨"



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com