



## VISCOMETER PETROLEUM EQUIPMENT BPTL-108

# VISCOMETER PETROLEUM EQUIPMENT BPTL-108

## AUTOMATIC CAPILLARY VISCOMETER WASHER

Viscosity of an oil is defined as its resistance to flow and Viscometer measures the flow. Monitoring viscosity is one of the most important components of any oil analysis program. Minute changes in viscosity at operating temperatures can be magnified to the extent that an oil is no longer able to provide adequate lubrication and efficiency in preventing friction between moving parts.



It adopts micro-controller control technology, program control, automatic washing. It needs less solvent oil and it is convenient, efficient. It can be widely used in the enterprises producing and using petroleum products, colleges and universities, and scientific research institutes and so on, where capillary viscometer are used. High degree automation.

## SPECIFICATIONS

Model	BPTL-108
Water bath temperature control range	Ambient to 80 °C
Ambient temperature	15 °C ~ 35 °C
Relative humidity	≤ 85 %
Maximum power consumption	700 W
Power supply	AC (220 ± 10%) V, 50 Hz
Dimension	320x300x500 mm (excluding the aspirator)
Net weight	24 kg (excluding the aspirator)
Alt Name	Automatic Capillary Viscometer Washer

## APPLICATIONS

Petroleum Industry, Oil and Gas Industry



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)