



## THERMO SHAKER INCUBATOR BSTH-206

# THERMO SHAKER INCUBATOR BSTH-206

## THERMO SHAKER INCUBATOR

Thermo shaker incubator designed for simultaneous heating, cooling and mixing of small samples. Both of them can be supplied with interchangeable blocks for microtubes.



Power failure self-recovery function, the device can automatically resume operation according to the originally set procedure after power failure recovery. Timing function, freely set the incubation time within 0~100 hours. Various standard sample blocks are available, which are easy to replace.

## SPECIFICATIONS

Model	BSTH-206
Temp. control range	RT. +5 °C~100 °C
Temp. uniformity	≤0.5 °C
Timer range	1 min~99 h 59 min
Temp. control accuracy	±0.5 °C
Display accuracy	±0.1 °C
Shaking speed	200-1500 rpm*
Heating time (40~100 °C)	≤10 min
Cooling time	/
Mixing orbit	2 mm
Dimension (mm)	300x220x170 mm
Weight	7 kg
Note	-
*The upper limit of temperature and shaking speed shall be determined according to the block used Laboratory Instrument	
Alt Name	Thermo Shaker Incubator

## ACCESSORIES FOR PURCHASE

No	Name	Description	Max. temp.	Max. speed
1	Block	96x0.2 mL tube	100 °C	1500 rpm
2	Block	54x0.5 mL tube	100 °C	1500 rpm
3	Block	35x1.5 mL tube	100 °C	1500 rpm
4	Block	35x2.0 mL tube	100 °C	1500 rpm
5	Block	15x0.5 mL+20x1.5 mL tube	100 °C	1500 rpm
6	Block	24xDia.≤Φ12 mm tube	100 °C	1200 rpm
7	Block	12x15 mL Falcon tube	100 °C	750 rpm
8	Block	6x50 mL Falcon tube	100 °C	750 rpm
9	Block	103x67x30 mm rectangular bath block	100 °C	1500 rpm
10	Block	96x0.2 mL microplate	100 °C	1500 rpm

11	Block	24x5 mL tube	100 °C	1200 rpm
12	Block	96 deepwell microplate	80 °C	800 rpm
13	Block	15x5 mL Eppendorf tube	100 °C	1200 rpm
14	Block	24x2 mL tube	100 °C	1200 rpm
15	Block	24x1.5 mL tube	100 °C	1200 rpm
16	Block	96x1000 µL deepwell plate	100 °C	1200 rpm
17	Block	MTPs deepwell plate	100 °C	1200 rpm



1



2



3



4




5



6



7



8



9



10



11



12



13



14



15



16



17



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)