



## STABILITY TEST CHAMBER BTST-1609

# STABILITY TEST CHAMBER BTST-1609

LARGE MEDICINE STABILITY TESTING CHAMBER



## SPECIFICATIONS

Model	BTST-1609
Temperature control range	0~65°C without Illumination 15~50°C with Illumination
Temperature Stability	±0.5°C
Temperature Uniformity	±2°C ~ ±3°C
Humidity Range	35~95%RH
Humidity Stability	±3%RH
Illumination	0~6000LX adjustable
Illumination difference	≤±500LX
Cooling system cooling mode	Two sets of Imported compressor work rotationally
Controller	Programmable (touch screen)
Sensor	Temp: Pt100, Humidity: capacitance sensor
Ambient Temperature	RT+5~30°C
Electrical Requirement	AC220V 50Hz
Power Consumption	3910W
Operating Power	1400W
Chamber Volume	500L
Interior Dimension (WxDxH,mm)	670x725x1020
External Dimension (WxDxH,mm)	850x1100x1930
Shelves	4(pcs)
Safety Device	Compressor overheating and overpressure protection, Fan overheating protection, Over temperature protection, Overload protection, Water protections
Remark	1. Build in mini printer 2. Paperless recorder (Option) 3. GSD series products have Installed Intensity of Illumination detector
Packing Dimension	124*97*213
CBM	2.56
Net Weight	317.00
Gross weight	400
Alt Name	Large Medicine Stability Testing Chamber



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)