



TUBE ROLLER / ROTATOR BTUR-401

TUBE ROLLER / ROTATOR BTUR-401

DIGITAL ROLLER MIXER

Digital Roller Mixer provides a gentle but highly efficient rocking and rolling motions, ideal for mixing blood samples, viscous substances and liquid-solid suspensions. Used in a variety of applications, prevention of blood coagulation and immune precipitation etc.



Rocking and rolling motions for complete mixing
Display shows the speed and time
Speed range of 10-99 rpm
Seven-roller design with small footprint and easy cleaning
Digital setting of speed and time

SPECIFICATIONS

Model	BTUR-401
Speed	Adjustable, 10-99 rpm ($\pm 10\%$)
Motion	Simultaneous roll & tilt ($\pm 5\%$)
Operating temperature	4-40 °C
Electrical	80-230 V, 50-60 Hz, 50 W
MOQ-10PCS	
Alt Name	Digital Roller Mixer

FEATURES

Rocking and rolling motions for complete mixing
Display shows the speed and time
Speed range of 10-99 rpm
Seven-roller design with small footprint and easy cleaning
Digital setting of speed and time



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com