



## HORIZONTAL AUTOCLAVE BAHZ-604

## HORIZONTAL AUTOCLAVE BAHZ-604

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## BAHZ-604 HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional

## SPECIFICATIONS

Model	BAHZ-604
Capacity	300 L
Chamber Size	Φ632x1000 mm
Chamber Material	SUS304
Max Working Pressure	0.23 Mpa
Working Temp.	115-138°C
Temperature Accuracy	+0.1° C
Consumption	17 KW
Power Supply	380 V, 50 HZ
External Size(LxWxH)	1400x890x1780 mm
Package Size(LxWxH)	1700x1050x2090 mm
Gross Weight	830 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)