



HORIZONTAL AUTOCLAVE BAHZ-801-A

HORIZONTAL AUTOCLAVE BAHZ-801-A

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

BAHZ-801-A HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional

SPECIFICATIONS

Model	BAHZ-801-A
Capacity	100 L
Chamber Size	Φ400x850 mm
Chamber Material	SUS304
Max Working Pressure	0.23 Mpa
Vacuum Limit	-0.08 Mpa
Working Temperature.	115-138°C
Temp Accuracy	10.1°C
Consumption	11 KW
Power Supply	380 V,50 Hz
External Size(Lx Wx H)	1250x730x1660 mm
Package Size(L x W x H)	1500x900x2050 mm
Gross Weight	310 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com