



## HORIZONTAL AUTOCLAVE BAHZ-801-B

## HORIZONTAL AUTOCLAVE BAHZ-801-B

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## BAHZ-801-B HORIZONTAL AUTOCLAVE



Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional

## SPECIFICATIONS

Model	BAHZ-801-B
Capacity	135 L
Chamber Size	Φ500x710 mm
Chamber Material	SUS304
Max Working Pressure	0.23 Mpa
Vacuum Limit	-0.08 Mpa
Working Temperature.	115-138°C
Temp Accuracy	10.1°C
Consumption	11.5 KW
Power Supply	380 V,50 Hz
External Size(Lx Wx H)	1110x750x1742 mm
Package Size(L x W x H)	1370x1000x2060 mm
Gross Weight	470 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)