



HORIZONTAL AUTOCLAVE BAHZ-801-E

HORIZONTAL AUTOCLAVE BAHZ-801-E

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

BAHZ-801-E HORIZONTAL AUTOCLAVE



- Sterilizing programs for fabric, instrument, rubber and liquid
- Touch screen display of temperature, pressure, time, running status and failure warning
- Post-drying function makes the residual humidity of sterilized articles less than 3%
- The positive pressure air-exhaust system to exhaust the cold air from chamber
- Equipped with printer to print date, time and process parameters
- Build-in efficient steam generator, microcomputer temperature control
- Perfect four-level authority management system to avoid misoperation
- Double door optional

SPECIFICATIONS

Model	BAHZ-801-E
Capacity	360 L
Chamber Size	Φ650x1130 mm
Chamber Material	S30408
Max Designed Pressure	-0.1/0.3 Mpa
Working Pressure	0.23 Mpa
Working Temperature.	105-134°C
Vaccum System	Water circulating vaccum pimp
Consumption	20 kVA
Noise	≤ 85 dB
Power Supply	380 V,50 Hz
External Size(Lx Wx H)	1570x1170x1980 mm
Package Size(L x W x H)	1755x1330x2270 mm
Gross Weight	770 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com