



HORIZONTAL AUTOCLAVE BAHZ-803-D

HORIZONTAL AUTOCLAVE BAHZ-803-D

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

BAHZ-803-D HORIZONTAL AUTOCLAVE

SPECIFICATIONS

Model	BAHZ-803-D
Capacity	241 L
Chamber Size (W x H x L)	600x600x740 mm
External Size (WxHxL)	1370x1880x1000 mm
Net Weight	720 Kg
Power Consumption	24 KW
Power Supply	380 V 50 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com