



## HORIZONTAL AUTOCLAVE BAHZ-804-F

## HORIZONTAL AUTOCLAVE BAHZ-804-F

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## BAHZ-804-F HORIZONTAL AUTOCLAVE

### SPECIFICATIONS

Model	BAHZ-804-F
Capacity	352 L
Chamber Size (W x H x L)	600x600x980 mm
External Size (WxHxL)	1215x1880x1190 mm
Net Weight	1000 Kg
Power Consumption	24 KW
Power Supply	380 V 50 Hz



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)