



## HORIZONTAL AUTOCLAVE BAHZ-810-B

## HORIZONTAL AUTOCLAVE BAHZ-810-B

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## BAHZ-810-B HORIZONTAL AUTOCLAVE



## SPECIFICATIONS

Model	BAHZ-810-B
Capacity	650 L
Chamber Size (W x H x L)	672x672x1500 mm
External Size (WxHxL)	1314x1260x1980 mm
Max Designed Pressure	Inner chamber - 0.1/0.3 Mpa, Jacket 0.3MPa
Working Pressure	0.23 Mpa
Built - in Steam Generator	0.8 Mpa
Designed Temperature	150 °C
Maximum Working Temperature	134 °C
Temperature Control Precision	± 0.2 °C
Consumption	56kVA
Standard Accessories	Trolley, Printer
Altitude	≤ 2000m; could be customised @ 2000-3000m
Power Supply	380 V 50 Hz
Package Size	Main Instrument: 1424x1460x2295; Trolley: 1860x640x1245
Gross Weight	Main Instrument: 1440; Trolley: 130



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)