



## CLEAN BOOTH BAIP-201

## CLEAN BOOTH BAIP-201

Clean booth is used as air purifier and it limits the contamination due to pollutants , thus providing partial high clean environment

Used in Pharmaceutical, Experimental drugs, Processing formula ,electronics, medicine, food, precision instruments, biochemical and chemical industries..

Also known as Down Flow booth.

## BAIP-201 CLEAN BOOTH



Customized design available

It can be used singly or combined

Compared with civil type and fabricated type clean room of hundred clean level, it has low running cost and fast effect and can be easily installed

Modular construction, easy to increase clean level, good expansion and reusable, convenient movement ( Universal wheel can be installed)

## SPECIFICATIONS

Model	BAIP-201
Clean Level	ISO 5(Class100),Class A
HEPA Filter	99.999% efficiency at 0.3µm
FFU Nos.	4 pcs
Control System	Microprocessor control system
Rectangular Tube Material	Aluminium alloy
Wall Material	Anti-static dustproof curtain
Blower	Built-in centrifugal blower, speed adjustable
UV Lamp	40W Emission of 253.7 nanometers with UV timer
Illumination Lamp	16Wx1 LED Lamp, Illumination: 2:500lux
Air Velocity	0.1-0.6m/s, average 0.45m/s
Internal Size (WxDxH)	1920x1380x2140 mm
External Size (WxDxH)	2000x1500x2530 mm
Package Size (WxDxH)	2110x930x920 mm 2280x660x500 mm
Consumption	800 W
Gross Weight	390 kg
Net Weight	355 kg
Noise	≤75 dB
Power Supply	AC220V±10%, 60Hz/50Hz; 110V±10%, 60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)