



Product Image Coming Soon

MEDICAL CLASS N AUTOCLAVE BAMN-404

MEDICAL CLASS N AUTOCLAVE BAMN-404

Designed for sterilization of solid instruments, medical class N autoclave is a non vacuum sterilizer which removes air from the chamber by thermodynamic air displacement. Safety features with door locking system, these autoclaves are ideal for small spaces.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

BAMN-404 MEDICAL CLASS N AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	BAMN-404
Capacity	23 L
Product Size	450x550x350 mm
Chamber Dimension	φ249x450 mm
Package Size	670x550x500 mm
Gross Weight	33 Kg
Net Weight	26 Kg
Power	800 W
Voltage	AC220/110V+10%, 50/60 Hz+ 1Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com