



HALOGEN MOISTURE ANALYZER BANA-503

HALOGEN MOISTURE ANALYZER BANA-503

Moisture Analyzer measures the amount of moisture in a substance which particularly helps in food processing. These analyzers are also commonly used in environmental, pharmaceutical, and materials science labs. We provide a variety of moisture analyzer according to applications, features and customers requirements.

Used in Food processing, Environmental, Pharmaceutical, Laboratory, Research, Medical, Hospitals, Agriculture, Industrial.

Also known as Laboratory Moisture Analyzer, Moisture Balance, Moisture Meter, Halogen Moisture analyzer, Laboratory Moisture Balance, Laboratory moisture Meter, Laboratory Halogen Moisture analyzer.

BANA-503 HALOGEN MOISTURE ANALYZER



- Halogen Light Heating
- Stainless Steel Chamber
- Temperature and time can be set
- The percentage of moisture content
- The percentage of dry residual
- The super high resolution back window matrix liquid crystal display (LCD), making it easier to operate in the dark place, and the user will have a more comfortable vision.
- Stores historical setting

SPECIFICATIONS

Model	BANA-503
Capacity	110g/1mg
Temperature	40°C-199°C
Operating Temperature Range	5°C-35°C
Temperature Set	1°C
Moisture Readability	0.0001
Dry Range	100.00%-0.00%
Dry Readability	0.0001
Pan Size	Φ90 mm
Readability	1 mg
Heat Source	HALOGEN LAMP-2
ATRO	100%...999%
ATRO 2	0%...999%
Sensor	LOADCELL
Temperature Sensor	PT-100
Time Setting	1-99 min By 10 s
Packaging Size	490x360x350 mm
Number of Storage	15
Weight	7.5 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com