



AIR SHOWER PASS BOX BAPB-101

AIR SHOWER PASS BOX BAPB-101

It is necessary passage for people entering clean room in-order to minimize particulate contamination
Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.
Also known as Pass Box.

BAPB-101 AIR SHOWER PASS BOX



Photoelectric sensor, automatic blowing
Adapt with circulation wind to ensure the cleanness of shower area under non-shower status
Double door by electric interlock
LED display with adjustable blowing time(0-99s)
HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BAPB-101
Clean Air Speed	0.36~0.54m/s
UV Lamp	15Wx1
External Size (WxDxH)	660x570x1000 mm
Internal Size (WxDxH)	500x500x500 mm
Material	304 stainless steel
Interlock	Electronical
HEPA Filter	Efficiency $\geq 99.995\%$ at 0.3 μ m
Gross Weight	80 kg
Packing Size(WxDxH)	780x700x1250 mm
Power Supply	220V, 50Hz; 110V, 50/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com