



AIR SHOWER PASS BOX BAPB-102

AIR SHOWER PASS BOX BAPB-102

It is necessary passage for people entering clean room in-order to minimize particulate contamination
Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.
Also known as Pass Box.

BAPB-102 AIR SHOWER PASS BOX

Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BAPB-102
Clean Air Speed	0.36~0.54m/s
UV Lamp	15Wx1
External Size (WxDxH)	760x670x1100 mm
Internal Size (WxDxH)	600x600x600 mm
Material	304 stainless steel
Interlock	Electronical
HEPA Filter	Efficiency ≥99.995% at 0.3µm
Gross Weight	100 kg
Packing Size(WxDxH)	880x800x1350mm
Power Supply	220V, 50Hz; 110V, 50/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com