



## PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-601

## PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-601

Reliable and robust class B autoclaves intended to deliver highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

## BAPV-601 PRE AND POST VACUUM CLASS B AUTOCLAVE

Class B, Pre & post vacuum

4 preprogrammed cycles, 1- customized programmed cycle

Integrated printer for tractability of sterilization cycles

Automatic regulation of the temperature

Light at the beginning of the cycle; Sound indication of end of cycle

No opening of the door during the cycle; Ability to interrupt the cycle at any time

Malfunction detection and warning system

Safety functions: auto-test program, over heat, over pressure, low water level indicator, safety door lock system present

Standard Accessory Sterilizing tray(2); Waste water tank(1); Tray removal tong(1); Tray holder(1)

## SPECIFICATIONS

Model	BAPV-601
Capacity	23 L
Sterilization Temperature	115°C, 121°C, 134°C
Working Pressure	-0.088 Mpa
Vacuum leakage rate	0.13 Kpa/min
Water quality	Distilled water
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Chamber Dimension	Ø250x450 mm
Overall Dimension (LxWxH)	620x470x385 mm
Display	LCD display
Weight	47 / 51 kg
Power	1900 W
Power Supply	220V/50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)