



PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-602

PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-602

Reliable and robust class B autoclaves intended to deliver highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

BAPV-602 PRE AND POST VACUUM CLASS B AUTOCLAVE

4 preprogrammed cycles, 1- customized programmed cycle

Integrated printer for tractability of sterilization cycles

USB port to save and download the recording of sterilization cycle

Built-in independent Steam generator to ensure the speed

Stand-alone drying system; Drying time can be edited

Safety functions: auto-test program, over heat, over pressure, low water level indicator, safety door lock system present

Standard Accessory Sterilizing tray(2); Waste water tank(1); Tray removal tong(1); Tray holder(1)

SPECIFICATIONS

Model	BAPV-602
Capacity	23 L
Sterilization Temperature	115°C, 121°C, 134°C
Water quality	Distilled water
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Fast sterilization cycle	16 min
Chamber Dimension	Ø250x450 mm
Overall Dimension (LxWxH)	620x470x385 mm
Display	5" color touch screen
Weight	48 / 52 kg
Power	1900 W
Power Supply	220V/50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com