



# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-702

## PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-702

Reliable and robust class B autoclaves intended to deliver highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

## BAPV-702 PRE AND POST VACUUM CLASS B AUTOCLAVE

Class B, Pre & post vacuum

4 preprogrammed cycles, 1- customized programmed cycle

Integrated printer for tractability of sterilization cycles

Automatic regulation of the temperature

Light at the beginning of the cycle; Sound indication of end of cycle

No opening of the door during the cycle; Ability to interrupt the cycle at any time

Malfunction detection and warning system, corresponding error code (F1-F19)

Double pressure release, suitable for sterilizing solid and liquid

Safety functions: auto-test program, over heat, over pressure, low water level indicator, safety door lock system present

Program 1 & 2: suitable for all instruments (wrapped or unwrapped, solid, hollow load products type A and porous products)

Program 3: suitable for solid, unwrapped products, single-layer wrapped products;

Program 4: Customized Program ( Sterilization Temperature, Vacuum time, Dry time and so on can be edited as required)

## SPECIFICATIONS

Model	BAPV-702
Capacity	45 L
Sterilization Temperature	115°C, 121°C, 134°C
Working Pressure	-0.098 Mpa
Vacuum leakage rate	0.13 Kpa/min
Water quality	Distilled water
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Chamber Dimension	Ø300x620 mm
Overall Dimension (LxWxH)	760x560x420 mm
Weight	86 kg
Power	3500 W
Power Supply	220V/50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)