



## PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-802

## PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-802

Reliable and robust class B autoclaves intended to deliver highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

## BAPV-802 PRE AND POST VACUUM CLASS B AUTOCLAVE

4 preprogrammed sterilization cycles, 1- customized programmed cycles  
Built-in Printer & USB port to save and download the recording of sterilization cycle  
Open-type water tank for easy cleaning  
Light at the beginning of the cycle; Sound indication of end of cycle  
No opening of the door during the cycle; Ability to interrupt the cycle at any time  
Malfunction detection and warning system, corresponding error code (F1-F19)  
Safety: Auto-test program, over heat/over pressure, low water level indicator, safety door lock system ; Safety for draining waste water

## SPECIFICATIONS

Model	BAPV-802
Capacity	60 L
Sterilization Temperature	115°C, 121°C, 134°C
Water quality	Distilled water
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Chamber Dimension	Ø320x700 mm
Overall Dimension (LxWxH)	950x600x580 mm
Display	5" color touch screen
Door	Automatic locking mechanism
Power	3800 W
Power Supply	220V/50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)