



## VERTICAL AUTOCLAVE BAVT-102

## VERTICAL AUTOCLAVE BAVT-102

We offer a variety of vertical laboratory autoclaves for general purpose sterilization applications. Customers can choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Food industry, Chemical industry, Research, Laboratory, Pharmaceutical.

Also known as Vertical Pressure Sterilizer, Vertical Pressure Autoclave, Vertical Sterilizer, Laboratory Vertical Autoclave,.

## BAVT-102 VERTICAL AUTOCLAVE



- Stainless steel structure
- Hand wheel type / Screw type
- Door safety lock system
- Digital Display
- Automatic exhaust system for cold air and pressure
- Over-temperature and over-pressure auto-protection
- Alarm function for completion of process and low water level
- Self-inflating type seal
- Optional: Drying system

## SPECIFICATIONS

Model	BAVT-102
Capacity	50 L
Temperature Range	105~134°C
Heat Average	±1°C
Working Pressure	0.22 MPa
Maximum Working Pressure	0.23 Mpa
Working Temperature	134 °C
Timer Range	0~90 min or 0~99hour59min
Overall Dimension (LxWxH)	520x490x1150 mm
Packing Size	580x580x1250 mm
Weight	94 / 74 kg
Power Supply	2x1.5KW/AC220V 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)