



PRECISION BALANCE BBAL-214

PRECISION BALANCE BBAL-214

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

BBAL-214 PRECISION BALANCE



Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing
4 conversion units: oz,ct,lb,g

SPECIFICATIONS

Model	BBAL-214
Capacity (g)	2100 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	≤2 sec
Pan Size	Φ 130 mm
Operation Temperature	5-25°C
Repeatability	0.02 g
Liner	0.03 g
Cal. Weight	1000 g
Packing Size	Inside: 305x215x115 mm Outside: 660x470x320 mm
Gross Weight	2 / 22 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer
3600715006	Windbreak mask



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com