



## PRECISION BALANCE BBAL-217

## PRECISION BALANCE BBAL-217

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

## BBAL-217 PRECISION BALANCE



Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing  
4 conversion units: oz,ct,lb,g

## SPECIFICATIONS

Model	BBAL-217
Capacity (g)	310 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	≤2 sec
Pan Size	Φ 130 mm
Operation Temperature	5-25°C
Repeatability	0.02 g
Liner	0.03 g
Cal. Weight	200 g
Packing Size	Inside: 310x230x130 mm Outside: 660x470x320 mm
Gross Weight	2 / 22 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer
3600715006	Windbreak mask



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)