



WEIGHING SCALE BBAL-417

WEIGHING SCALE BBAL-417

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory. Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research. Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

BBAL-417 WEIGHING SCALE



Stainless steel Pan
LCD digital display
The multiple functions modes with tare, counting, unit conversion and weighing

SPECIFICATIONS

Model	BBAL-417
Capacity (g)	61 kg
Resolution (g)	1 g
Min Weighing (g)	5 g
Stable Time	≤2 sec
Display	LED
Sensor	Load Cell
Pan Size	410x510 mm
Operation Temperature	5-25°C
Repeatability	±2 g
Liner	±3 g
Packing Size	420x170x850 mm
Gross Weight	13 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com