



## SHAKING WATER BATH BBSW-101

## SHAKING WATER BATH BBSW-101

Microprocessor based, digital shaking baths for gentle heating with agitation for your samples. Easy to clean stainless steel chamber, equipped with clamps that keep samples in place avoiding the risk of spills.

Used in Clinical, General Chemistry, Preservation, Reactions of sample, DNA amplification, Initial Denaturation of electrophoresis, Laboratory.

Also known as Digital Heating Block, Laboratory Mini Dry Bath Incubator.

## BBSW-101 SHAKING WATER BATH



Microprocessor based keypad for temperature/speed controls

Digital display

Easy to clean stainless steel internal chamber

Equipped with clamps to keep samples in place

It is featured with Spring Bottle rack, Alarm System, Adjustable Thermostat and Electronic Calibration

## SPECIFICATIONS

Model	BBSW-101
Capacity	50 L
Maximum Configuration	6 x 500ml or 3-4 x 1000ml
Rotation Range	0~220 rpm
Temperature Range	RT+5°C to 99°C
Time Range	0-999 h
Platform Dimension	380x300 mm
Outer Dimension (WxDxH)	705x515x470 mm
Weight	38 kg
Method	30 mm
Thermo Lid	Orbital
Power	1600 W
Power Supply	220V, 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)