



# PORTABLE LOW SPEED CENTRIFUGE BCBL-402

# PORTABLE LOW SPEED CENTRIFUGE BCBL-402

## BCBL-402 PORTABLE LOW SPEED CENTRIFUGE

Maintenance -free brushless motor  
LCD display of Speed, RCF and time  
Adjustable speed from 300r/min to 5000r/min  
Timer from 1s to 99min or continuous operation  
For 6x15ml centrifuge tube &10ml/7ml/5ml/3ml blood collection tube or 4x50ml centrifuge tube  
Lid automatically open at the end of run to revert sample warming and to reduce handling time  
10 Program stored in memory and editable by ender user  
10 Acceleration/ 10 Deceleration profiles ,with soft brake system to get the good separation effect  
Electric cover lock, compact design, super speed and imbalance protection

## SPECIFICATIONS

Model	BCBL-402
Maximum Capacity (No of tubes x Vol.)	4 x 50 ml
Maximum Speed	5000 rpm
Speed Accuracy	± 20 r / min
Maximum RCF	2635 × g
Time Range	0 ~ 99 min
Acc / Dec	10 Kinds
Noise Level	≤ 55 dBA
Certificates	CE,ISO & Calibration report are available
RPM / RCF Convert	Yes
Overall Dimension	280 x 220 x 175 mm
Net. Weight	3.5 kg
Temperature	Normal
Power Supply	AC 220 V /110 V 50 HZ / 60 HZ



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)