



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-101

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

BCBT-101

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

BCBT-101 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	BCBT-101
Capacity	4 gallon (15L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	430x430x560 mm
Package Size	550x560x725 mm
Gross Weight	38 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com