



## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-101

## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-101

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.  
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-101 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



## SPECIFICATIONS

|                  |  |
|------------------|--|
| Model            | BCBT-101   |
| Capacity         | 4 gallon (15L)                                   |
| Material         | Double layer anti fire cold-rolled steel PP tray |
| Adjustable Shelf | 1 pc   |
| Door             | Single / Manual                                  |
| External Size    | 430x430x560 mm                                   |
| Package Size     | 550x560x725 mm                                   |
| Gross Weight     | 38 Kg  |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)