



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-103

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

BCBT-103

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

BCBT-103 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	BCBT-103
Capacity	22 gallon (83L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	3 pcs
Door	Single / Manual
External Size	600x460x1650 mm
Package Size	720x580x1815 mm
Gross Weight	116 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com