



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-106

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-106

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.
Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

BCBT-106 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	BCBT-106
Capacity	60 gallon (227L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	860x860x1650 mm
Package Size	980x980x1815 mm
Gross Weight	195 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com