



## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-107

## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET BCBT-107

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-107 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



## SPECIFICATIONS

|                  |  |
|------------------|--|
| Model            | BCBT-107   |
| Capacity         | 90 gallon (340L)                                 |
| Material         | Double layer anti fire cold-rolled steel PP tray |
| Adjustable Shelf | 2 pcs  |
| Door             | Double / Manual                                  |
| External Size    | 1090x860x1650 mm                                 |
| Package Size     | 1210x980x1815 mm                                 |
| Gross Weight     | 243 Kg   |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)