



## STRONG ACID AND ALKALI STORAGE CABINET BCBT-202

## STRONG ACID AND ALKALI STORAGE CABINET BCBT-202

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-202 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti-corrosive materials for strong acid and alkali, Customizable

## SPECIFICATIONS

Model	BCBT-202
Capacity	12 gallon (45L)
Material	pp
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	590x460x890 mm
Package Size	675x545x1025 mm
Gross Weight	30 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)