



STRONG ACID AND ALKALI STORAGE CABINET BCBT-206

STRONG ACID AND ALKALI STORAGE CABINET BCBT-206

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

BCBT-206 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti-corrosive materials for strong acid and alkali, Customizable

SPECIFICATIONS

Model	BCBT-206
Capacity	45 gallon (170L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1650mm
Package Size	1175x545x1785 mm
Gross Weight	72 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com