



## STRONG ACID AND ALKALI STORAGE CABINET BCBT-207

## STRONG ACID AND ALKALI STORAGE CABINET BCBT-207

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## BCBT-207 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti-corrosive materials for strong acid and alkali, Customizable

## SPECIFICATIONS

Model	BCBT-207
Capacity	60 gallon (227L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	860x860x1650 mm
Package Size	945x945x1785 mm
Gross Weight	90 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)