



BLOCK CABINET BCBT-302

BLOCK CABINET BCBT-302

Three modules plus a base make up the cabinet. A combined stamping and shaping process is used to create labelling slots.

Used in Pathological laboratories, Research institutions, Universities, Cytology and Histology applications.

BCBT-302 BLOCK CABINET

The cabinet includes a base and three modules.

Labelling slots are manufactured through an integrated stamping and forming operation.

Hidden handles are made of chrome-coated ABS.

SPECIFICATIONS

Model	BCBT-302
Drawers	18 drawers (6 drawers in each module)
Capacity	11,000 tissue blocks
Base Material	SPCC1.2 cold-roll steel sheet
Main body and Drawers Material	SPCC 0.8 cold-roll steel sheets
Dimensions	450x478x1295 mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com