



SLIDE DRYING CABINET BCBT-303

SLIDE DRYING CABINET BCBT-303

Three modules plus a base make up the cabinet. The hidden handles are made of ABS. The slots are specifically manufactured using a metal mold.

Used in Pathological laboratories, Research institutions, Universities, Cytology and Histology applications.

BCBT-303 SLIDE DRYING CABINET

The cabinet includes a base and three modules.

Labelling slots are manufactured through an integrated stamping and forming operation.

Numbered labels inside the drawers
make filing more convenient.

Hidden handles are made of chrome-coated ABS.

SPECIFICATIONS

| Model | BCBT-303 |
|--------------------------------|---------------------------------------|
| Drawers | 18 drawers (6 drawers in each module) |
| Capacity | 9,000 microscope slides |
| Base Material | SPCC1.2 cold-roll steel sheets |
| Main body and Drawers Material | SPCC 0.8 cold-roll steel sheets |
| Dimensions | 450x478x1295 mm |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com