



AIR COOLED CHILLER BCHI-112

AIR COOLED CHILLER BCHI-112

Air cooled chiller absorb heat from process water, and the heat is then transferred to the air around the chiller unit. Air cooled chiller require less maintenance.

Used in Laboratory, Pharmaceutical, Chemical, Laser.

Also known as Laboratory Air Cooled Chiller.

BCHI-112 AIR COOLED CHILLER



It adopts single-stage vapor compression circuit and has compressor overload protection, pump overload protection, reverse phase and lack of phase warning, anti-icing protection, high and low pressure protection and other devices.

The machine has stable performance and long life.

It can cool down quickly, and the temperature is stable to meet customer requirements.

This series of products mainly work on the principle of cold and heat exchange.

It is suitable for the cooling field in modern industry and is not affected by the ambient temperature.

It is an indispensable configuration device.

SPECIFICATIONS

Model	BCHI-112
Freezing capacity	
kw	56.7 kw
Compressor	
Power	12.5hp×2
Condensor(Water Chiller)	
Type	Air cooling
Water tank capacity	250 L
Pump	
Power kw	3 hp
Flow rate	315 l/min
Weight	980 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com