



## AIR COOLED CHILLER BCHI-119

## AIR COOLED CHILLER BCHI-119

Air cooled chiller absorb heat from process water, and the heat is then transferred to the air around the chiller unit. Air cooled chiller require less maintenance.

Used in Laboratory, Pharmaceutical, Chemical, Laser.

Also known as Laboratory Air Cooled Chiller.

## BCHI-119 AIR COOLED CHILLER



It adopts single-stage vapor compression circuit and has compressor overload protection, pump overload protection, reverse phase and lack of phase warning, anti-icing protection, high and low pressure protection and other devices.

The machine has stable performance and long life.

It can cool down quickly, and the temperature is stable to meet customer requirements.

This series of products mainly work on the principle of cold and heat exchange.

It is suitable for the cooling field in modern industry and is not affected by the ambient temperature.

It is an indispensable configuration device.

## SPECIFICATIONS

Model	BCHI-119
Freezing capacity	
kw	128.6 kw
Compressor	
Power	12.5hp×4
Condensor(Water Chiller)	
Type	Air cooling
Water tank capacity	370 L
Pump	
Power kw	6 hp
Flow rate	400 l/min
Weight	1400 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)