



HIGH TEMPERATURE CIRCULATOR BCHT-104

HIGH TEMPERATURE CIRCULATOR BCHT-104

We offer you robust circulators for high temperature applications equipped with a circulating pump to provide excellent temperature uniformity and safe operation. Easy to maintain with LCD display, these circulators are user friendly and give reliable results.

Used in Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals.

Also known as Laboratory water bath.

BCHT-104 HIGH TEMPERATURE CIRCULATOR



Microprocessor temperature controller

Easy maintenance

LCD display with touch key panel for temperature setting

Over temperature alarm system

SPECIFICATIONS

Model	BCHT-104
Capacity	30 L
Temperature Range	RT+5°C to 100°C
Temperature Accuracy	±0.05°C
Pump Flow Rate	6 L/min
Top Opening	310x280 mm
Inner Dimension (WxDxH)	400x330x230 mm
Outer Dimension (WxDxH)	550x490x600 mm
Weight	21.62 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com