



MINI CENTRIFUGE BCMI-201

MINI CENTRIFUGE BCMI-201

Designed with quick spin option, our product is a compact, easy to use, high speed mini centrifuge that delivers maximum productivity. Variable speed allows you to carry out complex as well as simple protocols. High speed accuracy, low maintenance and safety operating features make it an ideal choice for your separation processes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.

Also known as Benchtop Centrifuge, Laboratory Centrifuge.

BCMI-201 MINI CENTRIFUGE



Brushless DC Motor ensures low maintenance, high efficiency, constant speed of rotation

Durable ABS housing resists chemical spills and limits noise

Short spin function is used for quick spin downs

Automatic release of lid when the rotor stops to save processing time

Safety Features: Over-speed detection, Over-temperature detection

SPECIFICATIONS

Model	BCMI-201
Maximum Capacity (No of tubes x Vol.)	6x1.5/0.5/0.2ml
Maximum Speed	4000 rpm
Speed Increment	±50 rpm
Maximum RCF	900xg
Time Range	None
Overall Dimension	170x150x130 mm
Display	None
Noise Level	≤ 48 dB(A)
Weight (Net/Gross)	1.25 kg
Power Supply	220 V/ 60 Hz

ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2302406006	Rotor	6x1.5/0.5/0.2ml	4000	900xg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com