



## FREEZE DRYER BFFT-313-A

## FREEZE DRYER BFFT-313-A

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption. Used in Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.

Also known as Floor Standing Freeze Dryer, Freeze Drying Equipment, Floor Type Freeze Dryer, Laboratory Floor Type Freeze Dryer, laboratory Floor Standing Freeze Dryer, Laboratory Freeze Drying Equipment, Floor Standing Lyophilizer, Floor Type Lyophilizer, T-Manifold Freeze Dryer.

## BFFT-313-A FREEZE DRYER



## SPECIFICATIONS

|                       |                                 |
|-----------------------|---------------------------------|
| Model                 | BFFT-313-A                      |
| Type                  | Top Press with Manifold Type    |
| Condenser Temperature | -110°C                          |
| Freeze Drying area    | : 0.09 m <sup>2</sup>           |
| Bulk capacity         | 0.9 L                           |
| Tray                  | φ200mm, 3pcs                    |
| Round bottom flask    | 1000ml,500ml ,250ml ,100ml,50ml |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)