



DUCTLESS FUME HOOD BFSD-504

DUCTLESS FUME HOOD BFSD-504

Ductless fume hoods by Biolab Scientific provide the needed protection from harmful gases, Vapours and toxic chemicals. Our ductless fume hoods have built-in airflows across the workstation, recycling clean air back into your laboratory. The Ductless fume hood is very effective if installed and used properly and maintained in good working order.

Used in Forensics, Life Sciences, Clinical, Research Field.

Also known as Carbon-filtered enclosures, Filtered fume hoods.

BFSD-504 DUCTLESS FUME HOOD



LCD Touch screen control panel, easy to operate.

With memory function in case of power-failure.

8°slope front ergonomics design, fatigue-free working posture.

Temperature and humidity sensors, can detect indoor temperature and humidity.

Three side transparent Acrylic windows, front window reversal design, easy to operate.

Double-layers structure: 1mm sheet metal surface; Chemical resistant phenolic resin work table.

Electronic control system, anti-overload, anti-electric shock, stable performance, long service life.

Inside and outside probe, detect indoor air pollution and filter conditions. Audible and visual alarm for changing filter.

SPECIFICATIONS

Model	BFSD-504
Air Velocity	0.4~0.6 m/s
Airflow Volume	547m³/h
Max.Opening	650mm
Work Surface Height	900mm
Blower	Built-in centrifugal blower, speed adjustable
Front Window	Acrylic window; Manual; Reversal design
Exterior Material	Cold-rolled steel with anti-bacteria powder coating
Side Window Material	Acrylic Window
Work Table Material	Chemical resistant phenolic resin
Chemical Filter	4 pcs
UV Lamp(optional)	Emission of 253.7 nanometers for most efficient decontamination
Fluorescent Lamp	14Wx3
Standard Accessories	Fluorescent Lamp, Base cabinet, Total load of 2 waterproof sockets: 500W;Chemical Filter
Optional Accessory	Water tap, Gas tap, Water sink, Electric height, Adjustable base stand, HEPA Filter, USB module,UV Lamp
Internal Size (WxDxH)	1710x690x740mm
External Size (WxDxH)	1800x940x2140mm
Main Body Package Size	1950x1110x1490mm
Package Size (Base Cabinet)	1950x1100x1060mm
Gross Weight	410kg
Consumption	500 W
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com