



-40°C UPRIGHT FREEZER BFUR-40-202

-40°C UPRIGHT FREEZER BFUR-40-202

Equipped with precise temperature controls, rugged interiors and safety features, our freezers provide you with the most dependable performance with low maintenance. A perfect partner for your multi cold storage purposes in a contemporary lab set up. Microprocessor based temperature controls provide outstanding temperature uniformity with high quality compressor for energy efficiency.

Used in Routine Cold Storage, Laboratories, Biologicals, Pharmaceuticals, Clinical, Research.
Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

BFUR-40-202 -40°C UPRIGHT FREEZER

- 4.3" color touch screen display
- Wide voltage design from 192V to 242V
- Pressure-balanced port design, easy to open the door
- Four walls and top with 90 mm, bottom with 130 mm
- Functions: High / Low Temperature, Hot condenser, Power Failure, Sensor Error, Low Battery, High Ambient temperature, Door Ajar, Remote Alarm
- Accessories: Caster, Test hole, Shelvesx3, Door locking device, Pressure balance hole, USB interface
- Optional Accessories: Chart Temperature Recorder, Racks & Boxes

SPECIFICATIONS

Model	BFUR-40-202
Capacity	480 L
Temperature Range	-10°C to -40°C
Cooling Performance	-40°C
Defrosting	Manual
Climate Class	N
Cooling Type	Direct Cooling
Refrigerant	Mixed Refrigerant, CFC free
Internal Dimension	600Wx580Dx1310H mm
Overall Dimension	780Wx780Dx1920H mm
Interior	Galvanized steel sheet



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com