



## -86°C ULT FREEZER UPRIGHT BFUR-86-317

## -86°C ULT FREEZER UPRIGHT BFUR-86-317

Freezers designed for daily sample protection and dependability. Environment friendly with safety features, these freezers provide outstanding thermal stability for a huge amount of samples. Modify your requirements with several optional features including chart recorder, CO2 back-up system, storage racks and remote alarm system.

Used in Preservation Of Plasma, Microbial Cultures, Tissues, Cell Lines, Clinical, Industrial, Research.

Also known as Laboratory ULT freezer, ULT Freezer, Ultra-low Freezer.

## BFUR-86-317 -86°C ULT FREEZER UPRIGHT



Cascade system, HC green refrigerant R290 & R170, Energy saving), Touch Screen, VIP

## SPECIFICATIONS

Model	BFUR-86-317
Capacity	728 L
Temperature Range	-40°C to -86°C
2ml Sample Loading	50400pcs
Alarm	High&low temp, power failure, control error, door open,abnormal voltage, low battery, condensor clean, high ambient temperature
Alarm port	USB+Remote alarm port
External Dimension	1229×1012×1994.5 mm
Optional Parts	110V transformer; Datalogger; Cold storage rack, box; LN2/CO2 backup system



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)