



AUTOMATIC AIR GENERATOR BGEN-102

AUTOMATIC AIR GENERATOR BGEN-102

Air Generator is designed for removing hydrocarbons from the air stream and produces pure, hydrocarbon free and dry air, with a very low dew point. Long service life, easy maintenance, abrasion resistant, cost effective, performance, quick and easy servicing are some of the features which makes it convenient for the laboratory.

Used in Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

Also known as Laboratory Air Generator.

BGEN-102 AUTOMATIC AIR GENERATOR



High precision output pressure

Water automatically drain out

Two levels stainless steel made purifications

Large inside stainless steel air tank

Have two levels pressure stable valve. Output air pressure much stable

Have over pressure protection, low voltage start up, DC, greatly protect inside compressor

SPECIFICATIONS

Model	BGEN-102
Output Pressure	0-0.4MPa
Noise	<35dB
Dimension	420x260x350 mm
Net Weight	22 kg
Consumption Power	150 W
Output Flow Rate	0-2000 ml/min
Environmental Temperature	0-45°C
Humidity	≤85%



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com